

## CARDIOVIT CS-200

Simply exceeds your expectations:

- Resting ECG
  - Vector Cardiography
  - Late Potential Analysis
  - QT Dispersion
- Exercise ECG
  - Ergospirometry
- Spirometry
- SEMA
- 24-h Holter ECG



# SCHILLER

The Art of Diagnostics

# Success through Innovation!

## CARDIOVIT CS-200



The CARDIOVIT CS-200 is a multitasking system solution, combining proven diagnostic capabilities with the modern computer technology and the most widely accepted user interface.

The CARDIOVIT CS-200 is the key to complete diagnostic services, including:

- 12-channel Resting ECG
- Automatic ECG measurements and interpretation (adult and pediatric)
- Vector Cardiography
- Pacemaker measurement
- QT Dispersion
- 12-channel full disclosure Exercise ECG, with ST monitoring and rhythm monitoring
- 24-h ECG Holter
- Fully integrated data management system
- Connection to network systems
- Resting rhythm recording
- Analysis of Ventricular Late Potentials
- Spirometry
- Ergospirometry

The open system architecture allows the implementation of new programs, future system expansion and network integration.

## Hardware Details



**DVD recorder: an easy way to update the system software or to archive the data.**



**Full-size keyboard with integrated trackball and diskette drive.**



**Direct function keys provide single button operation for immediate ECG printouts.**



**Integrated thermal printer or external printer of your choice.**

## Filter

- **SCHILLER** Powerline Filter (SPF):

Distortion-free suppression of superimposed 50 or 60 Hz sinusoidal interferences by adaptive digital filtering

- **SCHILLER** Smoothing Filter (SSF) / **SCHILLER**


**Myogram Filter (SMF-muscle tremor filter):**  
25 Hz/35 Hz low pass smoothing filter to suppress muscle tremor and high-frequency noise, without bothering the QRS complex

- **SCHILLER** Baseline Stabilizer (SBS):

Filter to suppress or greatly reduce baseline fluctuations without changing the measurement values (for Resting and Exercise ECGs)

# The innovative ECG System for the Diagnostic Professional!

## CS-200 Resting ECG



### Good Lead Quality saves Time

The exceptionally high digital sampling rate (2000 samples per second) recognizes every nuance of an ECG for maximum accuracy.

The CARDIOVIT CS-200 produces perfect ECG recordings every time by utilizing SCHILLER's unique filter technology and electrode hook-up test.

### Connect with the Future

The CARDIOVIT CS-200 can be integrated directly into a user network through standard technologies allowing for central archiving and data access from various locations.

### Ready to go at all Times


The integrated back-up battery guarantees basic system operation during power failures or in emergency situations. Automatic ECG recording and real-time printouts are available without delay.



The **CARDIOVIT CS-200** diagnostic workstation combines proven ECG diagnostic performance with the most advanced computer technology available.

- Built-in, full page thermal printer
- Real-time ECG recording
- Integrated data management system SEMA-200
- Windows™ user-interface
- Network capabilities

The unique modularity of the **CARDIOVIT CS-200** allows physicians to configure the system to meet their individual needs for professional diagnostics. In addition to resting ECG, the **CARDIOVIT CS-200** may be utilized to perform exercise testing and Holter monitoring with the highest level of user comfort and reliability.



The **CARDIOVIT CS-200**'s built-in thermal printer records 3, 6, or 12 channels simultaneously in real-time. The outstanding print resolution and wide selection of print formats make each report truly impressive.

Reports may even be generated on the optional laser printer using plain paper. Each **CARDIOVIT CS-200** comes with automatic measurements for all 12 standard leads.

The optional interpretation program generates automatic adult and pediatric interpretation statements. Interpretative statements can easily be edited before printing or archiving. Global measurements points can also be edited.

# Exercise Testing for the Professional!

## CS-200 Exercise Testing



View after exercise

Thermo printout

Averaged QRS complex

### Access the World of PC Technology

An electrocardiograph with the properties of a personal computer!

- Large color screen with XGA resolution for paperless operation
- Windows™ operating system
- Data management solution SEMA-200 for archiving and management of data
- Open system architecture for networking
- Common printer interfaces for documentation on standard (A4) paper
- Selection of built-in storage solutions (hard disk, diskette, CD-R/RW (option) etc.)
- Integration with existing management and information systems

### The CARDIOVIT CS-200 – the exemplary Exercise Testing Solution

The CARDIOVIT CS-200's exercise testing module fully demonstrates the system's capabilities. The large, color display provides clear monitoring throughout the entire test. Vital information at-a-glance for maximum patient safety and test reliability. Printouts are generated at pre-programmed intervals or at the simple push of a button.



The **CARDIOVIT CS-200** exercise testing system combines **SCHILLER's** proven ECG know-how with today's most modern, state-of-the-art personal computer innovations.

### Covering the full Spectrum of Electrocardiography

The **CARDIOVIT CS-200** offers the performance of a conventional exercise testing system:

- Built-in, full-size thermal printer for immediate real-time printouts
- No start-up delay – instant ECG reports at the push of a button
- Quick recording and printouts of ECGs
- Precise ECG measurements
- RS-232 and parallel interfaces
- Emergency ECG

**\*) Option: Blood Pressure Measurement BP-200 plus** is a non-invasive Stress Test Blood Pressure Monitor, designed to automatically measure and display a patient's systolic and diastolic pressures at preset intervals, or upon demand of the operator.

### The Convenience of automated Exercise Testing

The **CARDIOVIT CS-200** controls a selection of peripheral equipment (ergometers, NIBP devices and treadmills) through integrated interface protocols. Standard, and user defined exercise protocols are stored. These can be changed according to the situation and individual preferences as desired. Manual intervention to modify an ongoing protocol is possible at any time. From the beginning, the exercise test runs automatically so the user is able to fully concen-

trate on the patient without distraction.

### Continuous ST Analysis

The **EXEC** analysis program computes continuous ST measurement for all 12 leads.

### High Level of Recording Quality when it counts

Use of multiple digital filters, such as **SBS SCHILLER** Baseline Stabilizer or **SSF SCHILLER** Smoothing Filter, reduces muscle artifact


and baseline wander without changing the tracing fidelity.

### ECG Memory

All exercise tests are stored in the form of a detailed final protocol. The complete test is stored. The **CARDIOVIT CS-200** provides tools for posttest diagnostics, measurements and/or documentation of the recordings.

# Managing Information

## SEMA-200 Data Management



Patient record

Serial comparison

Search

### Simple to operate and learn

SEMA-200 is designed for daily use. Symbols, accessible by mouse, enable for quick learning and intuitive operation.

### Data Collection

SEMA-200 is installed directly onto the CARDIOVIT CS-200 workstation. Data can also be transmitted from any **SCHILLER**

product equipped with a communication interface.

### Data Management

Collected test data can easily be managed by patient name or ID, type of recording, including date and time or administrative status (validated, printed, sent, archived etc.).

### Serial Test Comparison

Multiple tests of the same patient can be compared side by side, an ideal tool for monitoring patient history or the effectiveness of therapy. Resting ECG only!

### Archiving

SEMA-200 stores acquired data forever and in any quantity. The program's compatibility with almost any current PC storage techno-



**SCHILLER's SEMA-200 data management program provides comprehensive management of diagnostic data from:**

- Resting ECG
- Exercise Testing
- QT Dispersion
- Spirometry
- Ergospirometry
- Analysis of Ventricular Late Potentials
- Vector Cardiography
- Holter ECG

**The program contains the following functions:**

- Validation with comprehensive editing tools
- Serial comparison
- Batch processing
- Archiving
- Network connection/integration
- Data communication (option with SEMACOM)
- Import/export for HIS integration

**SEMA-200 is a Windows™ XP Professional compatible application and may be run on the CARDIOVT CS-200 workstation, in a network or a stand-alone personal computer.**

logy makes it extremely flexible to fit a specific need. Storage on the internal hard disk or removable storage systems (optical disk, cartridge etc.) is possible.

### **Networking**

SEMA-200 fits seamlessly into networks providing access to data from different sites. Passwords and access limitations protect

the unauthorized manipulation of stored information. SEMA-200 runs on all current PC-based network solutions.

### **Data Transmission**

SEMA-200 communicates via direct cable or telephone modem to all **SCHILLER** equipment via the RS-232 interface. Data files can also be sent via the Internet.

### **Open System Architecture**


SEMA-200 is ready to be integrated into existing information systems.

# Estimating Arrhythmic Risk!

## CS-200 Late Potential Analysis

- Recognise Ventricular Late Potentials thanks to signal averaging and high pass filtering
- Gain time thanks to straightforward measurement data recording; a single key stroke suffices
- Comprehensive and clear evaluation within the shortest possible time
- QRS start and end of vector amplitude can be edited on-screen
- The optimum reference ECG lead for signal averaging can be defined automatically or can be selected by the user

Particularly for patients who have suffered a myocardial infarction, micro-potentials after the QRS complex indicate a higher risk of sudden cardiac death. The existence of Ventricular Late Potentials indicates an arrhythmic substrate that can cause re-entry circuits and may result in potentially lethal arrhythmia.



With the SCHILLER Late Potential analysis, it is possible to measure in Frank or in bipolar leads. The analysis is made within the time domain.

**Signal averaging\*** on a selectable quantity of heartbeats improves the signal/noise ratio. Together with a matured **high-pass filter technique\*** it enables recognition of the Ventricular Late Potentials.

The following **parameters\*** are determined automatically:


- High frequency QRS interval
- The RMS value of the last 40 msec.
- The portion of vector amplitudes below 40  $\mu$ V (LAS: low amplitude signal)

The limit values of these parameters identify whether or not a patient has Ventricular Late Potentials and are given in the acronyms of the interpretation.\*

\* The SCHILLER Late Potential analysis conforms to standards issued by the Task Force Committee of the European Society of Cardiology, the American Heart Association and the American College of Cardiology, in accordance with the publication entitled "Standards for Analysis of Ventricular Late Potentials Using High-Resolution or Signal-Averaged Electrocardiography", G. Breithard, M.E. Cain, N. El-Sheiff, N. Flowers, V. Hombach, M. Janse, M.B. Simson, G. Steinbeck. JACC 1991; Vol. 17 (5): 999-1006

# Non-invasive Method for Identifying the Risk of Lethal Arrhythmia!

## CS-200 QT Dispersion



Using the **CARDIOVIT CS-200** or other **SCHILLER** ECG machines and the **SEMA-200** diagnostic systems, it is possible to determine the spatial QT dispersion in the resting ECG. The software calculates the QT interval in each ECG lead.

The average of the QT intervals and the standard deviation of all QT intervals are calculated from the mean value.

In addition, the minimum and maximum QT intervals and the lead in which they occur are specified.

The program also calculates the difference of the QT interval in each lead to the calculated average QT interval ( $\Delta$ QT). And finally, the difference between the maximum and minimum QT interval is also analysed as  $\Delta$ max QT. This is considered as particular indicator for QT dispersion.

The QT times can be corrected manually for all leads.

QT dispersion indicates how inhomogeneous the ventricular repolarisation is. The greater the QT dispersion, the more inhomogeneous the ventricular repolarisation is and the greater the risk for a sudden cardiac death of the patient.

QT dispersion analysis can be used as a further non-invasive method – in addition to Late Potential analysis and heart rate

variability – for patients with an increased risk of life-threatening arrhythmia.

Recognise the non-invasive risk indicators of increased arrhythmic risk using the Late Potential analysis and QT dispersion analysis programs obtainable from **SCHILLER** – both can be run on the **CARDIOVIT CS-200**. Examine heart rate variability using the 24-hour **MT-101/200** ECG analysis system.


# The Innovative Diagnostic System!

## CS-200 Spirometry

In addition to:

- 12-channel Resting ECG
- 12-channel Exercise ECG
- 24-hour Holter ECG
- QT dispersion
- Vector Cardiography
- Late Potential analysis
- Integrated data management for all the above examinations
- The diagnostic station CARDIOVIT CS-200 now also offers pulmonary function testing as an option. This new possibility completes the clinical requirements of non-invasive diagnostic.

As always, **SCHILLER** uses the latest technology to simplify pulmonary function testing.



### Comprehensive Pulmonary Function Testing

Four measurement programs for inspiratory and expiratory tests are available: FVC, SVC, MVV and MV, with nominal/actual differences calculated automatically and pre/post medication comparisons to aid diagnostics.

### Easy Operation

One button push starts the tests. Tests can

be carried as many times as desired with the results clearly displayed on the large monitor. The flow volume trace is shown in real time and in color.

### Storage and Printing

Pulmonary function data is saved using the on-board data management SEMA-200; all previous examinations are easily accessed and viewed on the monitor.

### Interpretation/Diagnosis

The software gives diagnosis indications when requested. You can select firm terminology from a data bank (acronyms) and extend it with your own text.

### Long-Term Comparisons/Trends

When a series of recordings from one patient is available, trend traces can be shown.

## CS-200 Ergo-Spiro

### All in one:

- Simultaneous presentation of 12 ECG and spirometry flow curves with possibility to compare all parameters
- Real Breath-by-Breath gas determination of metabolic and hemodynamic parameters during rest and exercise tests
- Bi-directional flow sensor insensitive to moisture
- Two monitors for the testing system
- Limit monitoring of any parameter with alarm functions
- Automatic control of ergometer, using configurable profiles
- Compact system
- Economical, as all accessories are sterilizable and reusable
- Network ability of CARDIOVIT CS-200 Ergo-Spiro



### A Multitude of Functions – all in one Device!

The CARDIOVIT CS-200 Ergo-Spiro from SCHILLER: The complete function analysis of heart, lung, blood circulation and metabolism. With this all areas of cardiopulmonary diagnostics are covered. Whether for cardiology, pneumology, sports and industrial medicine, this multifunctional exercise testing system is the solution for every need.

### Breath-by-Breath Analysis!

Every breath and every heartbeat is analyzed according to the individually set parameters and shown on the monitors in real time. With its monitor the one-of-a-kind CARDIOVIT CS-200 Ergo-Spiro allows the direct visual comparison of ergo-spirometrical and ECG data! Various interfaces warrant perfect data exchange of the saved values to Patient Data Management Systems (PDMS).




# Template Matching in the MT-101/200 long-term ECG System

## MICROVIT MT-101/200 Holter System



Original size




The MT-200 program classifies automatically and with very high reliability 20 different ventricular and supraventricular events (couplets, triplets, bigeminy, trigeminus, absolute arrhythmia, tachycardia etc.). During template matching QRS complexes with the same shape are automatically grouped in a pattern template by the MT-200 program.

The program altogether forms, up to 12 "normal" templates, 1 "artefact" template, up to 300 different ventricular extrasystole or "VES" templates and up to 12 different supraventricular extrasystole or "SVES" templates, thereby providing you with a comprehensive summary of the heart's various stimulation types (spreads of stimulus).

With a click of the mouse you can examine the individual beats in the template. Re-classification of whole templates and individual QRS complexes is very simple: press the right mouse-button and select "Normal", "VES", "SVES", "Artefact" or "Delete" from the menu.



Drag



Draw



Drop



### Important:

You do not need previous PC experience in order to use the SCHILLER MT-200 analysis program.

### 'Drag&Drop' Function

You can combine one whole template with another: click on it, move to the appropriate target-template with the mouse-button depressed (drag) and release (drop). By this means the whole template is re-classified in two seconds; event analysis is modified automatically.

The **SCHILLER** MT-101/200 long-term ECG system has decisive advantages such as:

- Direct testing of ECG signal quality from the PC before recording starts
- The exceptionally small and light MT-101 Holter recorder
- WINDOWS™ NT/2000/xp-based MT-200 analysis software (state of the art)
- Summary showing trend in ST and heart-rate

The **SCHILLER** MT-101/200 long-term ECG system includes the following options:

- Analysis of 2 channels in different combinations
- ST Analysis
- Template Matching
- HR Trend overview
- HR Variability
- Pacemaker Templates
- Reclassification in less than 2 seconds with 'drag&drop' function

# CARDIOVIT CS-200

## Technical Data:

### System:

**Dimension:** 60 x 62 x 153 cm (l x w x h)

**Weight:** 71 kg

**Monitor:** Size: 17" TFT LCD

#### Operator interface:

- Standard alphanumeric PC keyboard with built-in trackball
- Mouse
- Direct function keys for Start/Stop/MAN/AUTO-ECG

**Power requirements:** 115/230 V (nominal), 50/60 Hz

**Power consumption:** Max. 80 VA

**Battery:** Built-in uninterruptible power supply (UPS) provides at least 3 minutes of back-up power to system and printer (without monitor) in the event of power failure.

**Operating system:** WINDOWS™ XP Professional multi lingual

#### Storage medium:

- Hard disk
- 3.5" diskette drive, 1.44 Mbytes
- DVD R/RW

### Hardware Options:

- Laser or inkjet printer
- ECG cable arm
- Electrode vacuum system
- Treadmill, ergometer or reclining ergometer
- BP-200 plus blood pressure device

All hardware options may be added at any time and are easily retrofitted.

### Printing:

**Chart paper:** Thermoreactive, z-fold, width: A4, approx. 60 m, ready-to-file

**Printing process:** Built-in high resolution thermal printer, 8 dots/mm – amplitude axis, 40 dots/mm – time axis, @ 25 mm/sec.

**Paper speed:** 5 / 10 / 12.5 / 25 / 50 mm/sec. (manual printout)

**Sensitivity:** 5 / 10 / 20 mm/mV, automatic or manual selection

**Print format:** 6/12 channel printout optimally shown on a width of A4 (200 mm), automatic baseline adjustment

**Channels:** 12

**External printer (optional):** Printout on plain paper with laser or inkjet printer

**Frequency response:** 0.05 – 150 Hz (IEC/AHA)

### Safety Standards:

**Specifications:** Patient input: Fully floating and isolated, defibrillation protected (only with original SCHILLER patient cable)

**Safety Standards:** IEC/EN 60601-1, UL 60601-1; C22.2 No. 601.1-M90, IEC/EN 60601-2-25, IEC/EN 60601-1-2 (EMC)

**Protection Class:** I according to IEC/EN 60601-1 (with internal power supply)

**Applied Part:** CF according IEC/EN 60601-1

**Conformity:** CE according Directive 93/42/EEC (Medical Devices)

**Classification:** IIa according Directive 93/92/EEC

#### Environmental conditions:

- Temperature, operating: 10° to 40° C
- Temperature, storage: -10° to 50° C
- Relative humidity, operating: 25 to 95 % (non-condensing)
- Pressure, operating: 700 to 1060 hPa

### Resting ECG:

- Simultaneous acquisition of all 8 active electrode signals from 12 leads
- Sampling frequency: 4000 Hz
- Pacemaker detection:  $\geq \pm 2$  mV/ $\geq 0.1$  msec.
- Emergency ECG
- Resting rhythm with event marking
- QT dispersion
- Review and re-measurement of all stored data

#### Options:

- Computer-aided ECG interpretation for pediatric and adult ECGs (C)
- Vector Cardiography
- Late Potential analysis (time domaine)

### Exercise ECG:

- Arrhythmia detection and review
- Standard and free programmable exercise protocols for bicycles and treadmills
- Automatic NIBP measurement
- Enlarged QRS complex with superimposed reference beat selectable for each lead
- Real time continuous ST analysis, amplitude and slope graphic
- Online ST measurement adjustment
- Re-analyse possibility of complete exercise test with ST, HR and NIBP trend graphic
- Borg scale, symptoms, end point criteria tables
- Emergency ECG any time possible

#### Options:

- Full disclosure rhythm review

### Filter:

**SPF SCHILLER Powerline Frequency Filter:** Distortion-free suppression of superimposed 50 or 60 Hz sinusoidal interferences by adaptive digital filtering

**SSF SCHILLER Smoothing Filter, SMF SCHILLER Myogram Filter (muscle tremor filter):** 25 Hz/35 Hz low pass smoothing filter to suppress muscle tremor and high-frequency noise, without bothering the QRS complex

**SBS SCHILLER Baseline Stabilizer:** Filter to suppress or greatly reduce baseline fluctuations without changing the measurement values (for Resting and Exercise ECGs)

### SEMA 200:

- Comprehensive management of diagnostic data of resting ECG, exercise ECG, spirometry, ergo-spirometry, 24-h ECG Holter, QT dispersion, LAPO
- Importing possibility via floppy disc, network or RS-232 interface
- Validation of interpretation, measurement values and patient data
- Serial comparison of resting ECG data and curves on the screen

### BP-200 plus Blood Pressure Measurement (Option):

#### System:

**Dimension:** 222 x 155 x 96 mm (l x w x h)

**Weight:** 1.08 kg (including SPO<sub>2</sub>, ECG Amp, excluding batteries)

**Power supply:** External supply+9 VDC output @ 0-1.7 amps

**LCD:** Backlit graphical LCD with multLanguage menu

#### Interface:

- ECG and RS-232 interface for SCHILLER ECG stress test systems
- USB interface

**Programming:** Menu guidance 5 buttons

**Safety standards:** IEC/EN 60601-1; UL 60601-1, C22.2

No. 601.1-M90; IEC/EN 60601-1-2 (EMC)

**Protection Class:** I according to IEC/EN 60601-1 (option: internal power supply)

**Applied Part:** BF/CF according IEC/EN 60601-1

**Conformity:** CE 0123 according Directive 93/42/ECC (Medical Devices)

**Classification:** II a according Directive 93/42/EEC

### Measurement:

#### Methods of measurement:

Primary: SCHILLER K-Sound Analysis with or without QRS trigger signal for stress test

Secondary: Oscillometric with smooth linear deflation for resting measurements

#### Measurement range:

- Systolic: 50 – 270 mmHg
- Diastolic: 20 – 150 mmHg
- Heart rate: 40 – 250 BPM

**Measurement intervals:** 2 – 20 min.

**Memory capacity:** 200 measurements

#### QRS Trigger:

- Signal from external ECG source

**Sampling intervals:** On command from ECG stress test system, manual or automatic measurement periods from 2 to 20 minute intervals

**Indicator for the user:** On screen message and audio tone for user defined ranges

### Optional unit version:

Built in SPO<sub>2</sub> module and internal ECG amplifier (to generate the QRS- trigger signal to support the blood pressure measurement)

### Standard accessories:

- Cuff (M)
- User guide
- Velcro fixing pads to secure cuff
- Adhesive pads to fix the microphone
- Ext. power supply
- Head phone
- Communication cable to stress test system
- QRS-trigger cable

### Optional accessories:

- Several cuff sizes
- Batterie charger
- Rechargeable batteries
- SPO<sub>2</sub> with cable and ear-sensor
- Pat. cable for internal ECG Amplifier (IEC coding, Code 1 / AHA coding, Code 2)
- Mini SD card for SW-Update and additional storage of the raw data

# Technical Data:

## Spirometry:

### Measured values:

**FVC:** FVC, FEV<sub>0.5</sub>, FEV<sub>1.0</sub>, FEV<sub>3.0</sub>, FEV<sub>0.5</sub>/FVC, FEV<sub>1.0</sub>/FVC, FEV<sub>3.0</sub>/FVC, FEF<sub>0.2-1.2</sub>, FEF<sub>25-75%</sub>, FEF<sub>75-85%</sub>, PEF, MEF<sub>75%</sub>, MEF<sub>50%</sub>, MEF<sub>25%</sub>, FIVC, FIV<sub>1.0</sub>, FIV<sub>1.0</sub>/FVC, FIV<sub>1.0</sub>/FVC, PIF, MIF<sub>50%</sub>

**SVC:** SVC, ERV, IRV, TV

**MV:** MV, RR, TV

**MVV:** MVV, RR, TV

### Presentation:

- Flow/volume loop
- Volume/time graph
- Measurement values table
- Trend display at >1 recording per patient
- Interpretation quadrant

### Prediction equation:

**Adults:** ECCS, Austria, Crapo, Morris, Knudson, Knudson 76, Polgar, Berglund, Finland, India, Composite

**Children:** Quanier & Tammeling, Austria, India, Knudson, Knudson 76, Polgar

### Extrapolated predicted values

**Comparison:** Pre/post medication possible

**Standards compliance:** ATS, ATS 94, OSHA, NIOSH

### SPIROVIT SP-150 Pneumotach flowsensor for pulmonary function testing with disposable mouthpiece:

**Dimensions:** 118 x 36 x 28 mm; approx. 120 g

**Measuring method:** Pneumotachometer

**Measuring accuracy:** According to ATS standards < 3 %

**Flow impedance:** < 0.2 mbar\* s/l at 12 l/s

### SPIROVIT SP-20 Pneumotach flowsensor for pulmonary function testing with reusable mouthpiece:

**Dimensions:** 125 x 36 x 28 mm, approx. 160 g

**Measuring method, measuring ranges, measuring accuracy, flow impedance same as SP-150.**

### Basic unit:

- CS-200 Spirometry program with 4 PC based measurement programs for inspiratory and expiratory pulmonary function tests
- 1 pneumotach sensor (choice of SP-150 for disposable or SP-20 for reusable mouthpieces)
- Accessories: 2 noseclips and either 1 pack of disposable plastic mouthpieces (for SP-150) or 1 pack of disposable filters (for SP-20)
- 1 operating manual

**Hardware option:** Calibration syringe, P/N 2.100027

## Ergo-Spirometry:

Technical Description Gas Analyzer:

**Dimensions:** (integrated into CS-200): 15 x 15 x 15 cm (l x h x w)

**Weight:** 2.35 kg

**Power supply:** 100 – 240 V, 50/60 Hz

**Safety Class:** VDE II, safety insulated, Class BF

Measurement Methods:

### Breath-by-breath gas analysis

#### Flow volume processing according to ECCS or ATS

#### Amount of gas adjusted with STPD

#### Environmental sensors:

- Flow sensor insensitive to moisture: 0 to ± 17 l/sec., resolution 10 ml
- Temperature: 0 to 50° C
- Barometer pressure: 600 to 1050 mBar (measurement range: -100 to 3500 m above sea level)

#### Flow measurement (patented diaphragm spirometer):

Bi-directional pneumotach flow sensor insensitive to moisture  
Measurement range 0 to ± 17 l/sec., max. deviation ± 2.5 %, volume range 0 to 300 l, resolution 10 ml/sec.

#### Gas analyzer

##### O<sub>2</sub> chemical cell:

Measurement range: 0 to 100 % O<sub>2</sub>  
Measurement accuracy: Better than ± 0.1 %  
Response time: < 80 msec. (10 to 90 %)

##### Ultra sound density method analyzer CO<sub>2</sub>:

Measurement range: 0 to 17.5 % CO<sub>2</sub>  
Measurement accuracy: Better than ± 0.1 %  
Response time: < 80 msec. (10 to 90 %)

#### Fully automatic zero alignment and calibration of the gas analyzers at minimal gas consumption

#### Automatic measurement of ambient conditions

#### Manual entry of blood gas values (lactate etc.)

Software:

**Parameters:** METS, VT, VE, fR, VO<sub>2</sub>, VCO<sub>2</sub>, RQ, VE/V, VO<sub>2</sub>, VCO<sub>2</sub>-equivalents, O<sub>2</sub>-pulse, VD/VT, FEO<sub>2</sub>, FECO<sub>2</sub> as well as secondary parameters

#### Automatic determination of the anaerobic threshold after V-slope, CO<sub>2</sub> excess, EQO<sub>2</sub> minimum or RQ

#### Scaling optimization for the diagrams

#### Editing possibility for all basic parameters

#### Export/import function with standardized protocols for integration into HIS

Operating Conditions:

**Ambient temperature:** Between 10 and 40° C

**Relative humidity:** Between 30 and 90 %

**Atmospheric pressure:** Between 900 and 1060 hPa

The PowerCube produced by Ganshorn Medizin Electronic GmbH, Niederlauer, Germany, has been integrated into the CARDIOVIT CS-200 Ergo-Spiro device. It fully complies with SCHILLER's view on quality.

## 24-h ECG Holter MT-101/200:

ECG acquisition:

### ECG amplifier:

- Galvanically isolated: 5.5 kVDC
- Input impedance: > 10 MΩ
- Maximum electrode potential: ± 300 mVDC
- Common mode rejection: > 80 dB
- Dynamic range: ± 5.12 mVAC
- Frequency response: 0.05 – 150 Hz (-3 dB)
- Digital resolution: 2.5 μV, 12 bit
- Sampling frequency: 500 or 1000 Hz
- Time constant: = 3.2 sec.
- Pacemaker detection: ± 2.200 mV / 0.1..2 msec.

### Patient cable:

- 2-lead: 4 electrodes
- 3-lead: 6 electrodes

Recorder:

**Dimensions:** 94 x 61 x 20 mm, approx. 110 g (including battery)

**Data recording:** Simultaneous 2- or 3-channel recording, differential

**Power consumption:** ~ 95 mW (recording)

### Battery:

- 1 pc. AA alkaline 1.5 V (Energizer or equivalent)
- Rechargeable accumulator, 1 x AA NiMH with 1.2 V and capacity ≥ 2000 mAh

**Storage medium:** 256 MB SD Memory Card

**Data transmission from MT-101:** With integrated USB interface or removable SD card

### Software features MT-200:

- Save in PDF-Format (Portable Document Format)
- E-Mail-Function
- Multiple analysis
- Connection to the SCHILLER SEMA-200 Data Management System
- HR trend
- ST analysis (option)
- HR variability (option)
- Template matching (option)
- Pacemaker templates (option)
- P-wave detection (option)
- Artefact detection (option)

All technical data is subject to changes aiming to continuous innovations.

**Asia**

**SCHILLER Asia-Pacific**  
7 Jalan Ara SD 7/3B  
Bandar Sri Damansara,  
52200 Kuala Lumpur  
Malaysia  
Phone +603 6272 3033  
Fax +603 6272 2030  
sales@schiller.com.my  
www.schiller-asia.com

**Austria**

**SCHILLER Handelsgesellschaft mbH**  
Kampmüllerweg 24  
A-4040 Linz  
Phone +43 732 70 99-0  
Fax +43 732 757 000  
sales@schiller.at  
www.schiller.at

**C.I.S. & Russia  
SCHILLER AG**

Moscow Representative Office  
Office 414, 1-st Yamskogo Poly Street 15  
125124 Moscow, Russia  
Phone/Fax  
+7 (495) 970 11 33,  
+7 (495) 956 29 10  
mail@schiller-ag.com  
www.schiller-cis.com

**France (office Paris)  
SCHILLER Médical S.A.S.**

7, rue Raul Follereau  
Bussy Saint Georges  
F-77608 Marne la Vallée/Cedex  
Phone +33 1 64 66 50 00  
Fax +33 1 64 66 50 10  
info@schiller.fr  
www.schiller-france.com

**France  
SCHILLER Médical S.A.S.**

4, rue Louis Pasteur  
F-67162 Wissembourg/Cedex  
Phone +33 3 88 63 36 00  
Fax +33 3 88 94 12 82  
info@schiller.fr  
www.schiller-france.com

**India**

**SCHILLER Healthcare India Pvt. Ltd.**  
D.C. Silk Mills Compound, 'A' Wing,  
1st floor  
5, Chunawala Estate, Kondivita Lane  
Andheri-Kurla Road, Andheri (E)  
Mumbai - 400 059, India  
Phone +91 22 2826 3520  
Fax +91 22 2826 3525  
mail@schillerindia.com  
www.schillerindia.com

**Japan**

**SCHILLER Japan, Ltd.**  
Biomedical Research Building 218  
1-2-3 Kasumi, Minami-ku  
Hiroshima 734-8551  
Phone +81 82 250 2055  
Fax +81 82 253 1713  
koji.maekawa@schiller.jp  
www.schiller.jp

**Spain**

**SCHILLER ESPAÑA, S.A.**  
Ctra. de la Coruña  
A-6 km.18,200  
E-28230 Las Rozas/Madrid  
Phone +34 91 713 01 76  
Fax +34 91 355 79 33  
schiller@schiller.es  
www.schiller-spain.com

**Switzerland**

**SCHILLER-Reomed AG**  
Riedstrasse 14  
CH-8953 Dietikon  
Phone +41 44 744 30 00  
Fax +41 44 740 37 10  
sales@schiller-reomed.ch  
www.schiller-reomed.ch

**Turkey**

**SCHILLER Türkiye**  
Halil Rifat Pasa Mah.  
Darülaceze Cadd.  
Perpa Ticaret Merkezi  
B/Blok No: 41 K: 2  
Okmeydanı-Sisli - İstanbul  
Phone +90 212 210 8681 (pbx)  
Fax +90 212 210 8684  
sales@schiller-turkiye.com  
www.schiller-turkiye.com

**Germany**

**SCHILLER Medizintechnik GmbH**  
Rudolf-Diesel-Strasse 14  
D-85521 Ottobrunn  
Phone +49 89 62 99 81-0  
Fax +49 89 609 50 90  
info@schillermed.de  
www.schillermed.de

**USA**

**SCHILLER America Inc.**  
2131 N.W. 79th Avenue  
Doral, Florida 33122  
Phone +1 786 845 06 20  
Fax +1 786 845 06 02  
sales@schilleramerica.com  
www.schilleramerica.com



# SCHILLER

The Art of Diagnostics