

cardiofax GEM

Interpretive Electrocardiographs

ECG-9010K
ECG-9020K

Everything you've always wanted in an ECG



ECG-9020K

VALUE • RELIABILITY • CUSTOMER SERVICE



Recording Example

12 Lead ECG

4-channel as shown

6-channel selectable



ECG-9010K

Model ECG-9010k
excludes display

OPTIONS

- NK-POUCH90 Carry Case
- Pc CARDIO Review

This brochure may be revised or replaced by Nihon Kohden at any time without notice.

This device is intended for use only by qualified medical personnel.



NIHON KOHDEN AMERICA, INC.
90 Icon Street
Foothill Ranch, CA 92610
Phone: (800) 325-0283/(949) 580-1555
Fax: (949) 580-1550

Document Control # 90112-A

Printed in U.S.

A Tremendous Value in a Small Package

Space-saving design

The Cardiofax GEM is a convenient answer for offices with limited space. The compact design requires minimal table surface for use and storage.

Packed full of value


The GEM offers a full range of capabilities for an affordable price, by incorporating a high-level analyses program, 12-lead and rhythm assessment, rechargeable battery power operation and optional 12-lead preview display.

Rapid assessment

The ECAPS 12C assesses waveforms promptly for high-level evaluation and measurements. This package will minimize countless hours of input and gives you the clear advantage of quick feed back.


Quality recordings

An integrated display (MODEL ECG-9020k) allows reviewing of all 12-leads prior to recording, ensures quality results.



Data Transfer

The GEM offers the ability to transfer ECG recording to a computer for physician remote review and archiving. A single channel waveform may be imported or exported from either model.



Selection Table

	ECG-9010K	ECG-9020K
12-Lead Recording	x	x
3,4,6 Recording Formats	x	x
Rhythm Page	x	x
Periodic Recording	x	x
Interpretive Analyses	x	x
Waveform Preview		x
Computer Storage (clinical)		x
AC and Battery Power	x	x

ECG-9020K (3/4 actual size)